

Account #: _____ Ship To #: _____ Date: _____

P.O. #: _____ Contact Name: _____ Phone #: _____
Patient First Name: _____ Patient Last Name: _____

☐ Use previous CAMISHA impression scan(s) to manufacture this order. Right Serial #: _____ Left Serial #: _____

☐ UNIQUE440 ☐ UNIQUE330 ☐ UNIQUE220 ☐ UNIQUE110

ZEN Program comes standard in all UNIQUE models.

CIC MICRO (RC-Coil and T-Coil not included): ☐ R ☐ L **NOTE: Requires WidexLink™ Connector for programming.**

CIC (RC-Coil included, T-Coil not included): ☐ R ☐ L

Custom (XP) Instruments (T-Coil and RC-Coil standard; VC and Program Button **not** available):

Canal ☐ R ☐ L Half Shell ☐ R ☐ L Full Shell ☐ R ☐ L

Standard CIC: CAMISHA Laser Fit Technology (hypoallergenic), Medium Brown Hearing Instrument, NanoCare™ Removal Line.

Standard ITE: CAMISHA Laser Fit Technology (hypoallergenic), Beige Hearing Instrument, NanoCare™.

NON-STANDARD OPTIONS**Hearing Instrument Colors**

☐ Beige ☐ Medium Brown ☐ Dark Brown

☐ Clear ☐ Faceplate Colors (If a clear shell is chosen, select a faceplate color)
☐ Clear ☐ Beige ☐ Medium Brown ☐ Dark Brown

Battery Door Colors (If a standard clear faceplate is chosen, select battery door color)

☐ Beige ☐ Medium Brown ☐ Dark Brown

Additional Options (Where available)

☐ Soft Hypoallergenic Shell Coat ☐ Removal Line ☐ Extended Canal Lock
☐ Hard Hypoallergenic Shell Coat ☐ Removal Notch (N/A - CIC) (N/A - Full Shell/CIC)
☐ Retention Ring ☐ Safety Loop (N/A - CIC) ☐ No Helix - 3/4 Shell
☐ Removal Post (N/A - CIC) ☐ Canal Lock (N/A - Full Shell)

Venting Options

If no vent option is selected, Widex will select the appropriate vent based on the 500 Hz Audiometric Data provided.

☐ Small (1.0 mm) ☐ IROS Vent ☐ Pressure Vent
☐ Medium (1.5 mm) ☐ Semi IROS Vent ☐ Largest Vent Possible
☐ Large (2.0 mm) ☐ Reverse Horn Vent ☐ No Vent
☐ Trench Vent

WARRANTY OPTIONS (See UNIQUE Price Sheet for price information.)

INCLUDED WARRANTIES	OPTIONAL WARRANTIES
UNIQUE440 3-Year Repair Warranty 3-Year Loss & Damage Coverage (Standard)	
UNIQUE330 3-Year Repair Warranty 2-Year Loss & Damage Coverage (Standard)	3 rd Year of Loss & Damage Coverage <input type="radio"/> R <input type="radio"/> L
UNIQUE220 2-Year Repair Warranty 1-Year Loss & Damage Coverage (Standard)	3 rd Year Repair Warranty <input type="radio"/> R <input type="radio"/> L 2 nd Year of Loss & Damage Coverage <input type="radio"/> R <input type="radio"/> L 3 rd Year of Loss & Damage Coverage <input type="radio"/> R <input type="radio"/> L
UNIQUE110 1-Year Repair Warranty 1-Year Loss & Damage Coverage (Standard)	2 nd Year Repair Warranty <input type="radio"/> R <input type="radio"/> L 3 rd Year Repair Warranty <input type="radio"/> R <input type="radio"/> L 2 nd Year of Loss & Damage Coverage <input type="radio"/> R <input type="radio"/> L 3 rd Year of Loss & Damage Coverage <input type="radio"/> R <input type="radio"/> L

DEX Items**Quantity**

COM-DEX (Specify color: Champagne White, Anthracite Grey, Emerald Green)

COM-DEX Remote Mic (Grey)

CALL-DEX (Black)

UNI-DEX (Black)

PHONE-DEX (Silver)

RC-DEX (Specify color: Silver, Black)

TV-DEX (Silver)

TV-DEX (Additional Base - Silver)

FM+DEX (Silver)

T-DEX (Silver) Not compatible with CIC

NOTE: CIC MICRO is not compatible with any DEX device.

Canal Length Options

☐ Medium (Standard)

☐ Short

☐ Long

☐ As Marked on Impression

Wax Guard Options

☐ NanoCare™ (Standard)

☐ Extended Receiver Tubing

☐ Spring

☐ Bell Canal

Impression Recommendations:

Silicone impression material; include second bend.

Shipping Method & 48-Hour Service

(See Price & Policy Guide for price information.)

☐ FedEx Next Afternoon (Standard)

☐ FedEx Delivery (Extra charge)

☐ 48-Hour Rush Service (Extra charge)

Please Send:

☐ Order Forms ☐ Impression Boxes ☐ Service Forms

AUDIOMETRIC DATA

	250Hz	500Hz (Required)	1kHz	2kHz	4kHz
Right					
Left					

Previous User ☒ Y ☐ N Age: _____

Ear Texture ☐ Soft ☐ Medium ☐ Firm

Additional Comments / Special Instructions: _____

