

WIDEX UNIQUE IM AND IP ORDER FORM



CUSTOMER INFORMATION

Bill-to Account Number: _____
 Ship-to Account Number: _____
 PO #: _____
 Address: _____
 City: _____ State: _____ Zip: _____
 Contact Name: _____
 Phone: _____ Email: _____

IMPRESSION RECOMMENDATIONS



- Use silicone impression material.
- Make impressions past the second bend.
- For optimal microphone placement, it is recommended to draw a horizontal line through the tragus on the earmold impression.

PATIENT INFORMATION

Date: _____
 First Name: _____
 Last Name: _____

Audiometric Information

If no vent option is selected on this order form, Widex will select the appropriate vent based on the Audiometric Data provided.

	250Hz	500Hz (Required)	1kHz	2kHz	3kHz	4kHz
Right						
Left						

Widex may decrease the vent size to accommodate your order request. Would you like to be consulted before this change is made? ☐ Y ☐ N

Widex may change your selected sides for the volume control and program button features. Would you like to be consulted before this change is made? ☐ Y ☐ N

Previous User: ☐ Y ☐ N Age: _____
 Ear Texture: ☐ Soft ☐ Medium ☐ Firm
 Impression: ☐ Open mouth ☐ Closed mouth

☐ Use previous CAMISHA impression scan(s) to manufacture this order. Right Serial #: _____ Left Serial #: _____

SELECT CUSTOM ITE OPTIONS

Model	Design	Hearing Instrument Colors	Venting May vary due to ear canal size	Additional Options	Canal Length
☐ UNIQUE440 ☐ UNIQUE330 ☐ UNIQUE220 ☐ UNIQUE110 ZEN Program comes standard in all UNIQUE models. (RC coil included. T-coil not included. 312 battery standard.)	Canal ☐ R ☐ L Half Shell ☐ R ☐ L Full Shell ☐ R ☐ L	☐ Beige ☐ Medium Brown ☐ Dark Brown ☐ Clear (If Clear shell is chosen, select a faceplate color) ☐ Beige ☐ Medium Brown ☐ Dark Brown	IROS Vent ☐ R ☐ L Semi IROS Vent ☐ R ☐ L Reverse Horn Vent ☐ R ☐ L Trench Vent ☐ R ☐ L Pressure Vent ☐ R ☐ L Largest Vent Possible ☐ R ☐ L No Vent ☐ R ☐ L XS ☐ R ☐ L S ☐ R ☐ L M ☐ R ☐ L L ☐ R ☐ L XL ☐ R ☐ L	Soft Hypoallergenic Shell Coat ☐ R ☐ L Hard Hypoallergenic Shell Coat (Default) ☐ R ☐ L Retention Ring ☐ R ☐ L Removal Line ☐ R ☐ L Removal Notch ☐ R ☐ L Canal Lock (N/A - Full and Half Shell) ☐ R ☐ L No Helix - 3/4 Shell ☐ R ☐ L	Medium (Standard) ☐ R ☐ L Short ☐ R ☐ L Long ☐ R ☐ L As Marked on Impression ☐ R ☐ L Wax Guard NanoCare™ (Standard) ☐ R ☐ L Extended Receiver Tubing ☐ R ☐ L Spring ☐ R ☐ L
Receiver	Feature*				
U-IM (85 dB) ☐ R ☐ L U-IP (90 dB) ☐ R ☐ L	Volume Control ☐ R ☐ L Program Button ☐ R ☐ L *To optimize size and ease of use, Widex recommends NOT adding a push button and volume control on both ears.				

WARRANTY OPTIONS (See UNIQUE Price Sheet for price information.)

INCLUDED WARRANTIES	OPTIONAL WARRANTIES
UNIQUE440 3-Year Repair Warranty 3-Year Loss & Damage Coverage (Standard)	
UNIQUE330 3-Year Repair Warranty 2-Year Loss & Damage Coverage (Standard)	3 rd Year of Loss & Damage Coverage ☐ R ☐ L
UNIQUE220 2-Year Repair Warranty 1-Year Loss & Damage Coverage (Standard)	3 rd Year Repair Warranty ☐ R ☐ L 2 nd Year of Loss & Damage Coverage ☐ R ☐ L 3 rd Year of Loss & Damage Coverage ☐ R ☐ L
UNIQUE110 1-Year Repair Warranty 1-Year Loss & Damage Coverage (Standard)	2 nd Year Repair Warranty ☐ R ☐ L 3 rd Year Repair Warranty ☐ R ☐ L 2 nd Year of Loss & Damage Coverage ☐ R ☐ L 3 rd Year of Loss & Damage Coverage ☐ R ☐ L

SELECT ACCESSORIES (Indicate quantity)*

RC2-DEX: Silver _____ Black _____ COM-DEX: White _____ Grey _____ Green _____ COM-DEX Remote Mic: Grey _____	UNI-DEX: Black _____ FM+DEX: Silver _____ TV-DEX: Silver _____ ADD'L TV-DEX BASE: Silver _____ PHONE-DEX: Silver _____ CALL-DEX: Black _____
--	---

*All DEX devices are billable items for IM/IP models

Additional Comments / Special Instructions: _____