

## CUSTOM IM AND IP ORDER FORM



### CUSTOMER INFORMATION

Bill-to Account Number: \_\_\_\_\_  
 Ship-to Account Number: \_\_\_\_\_  
 PO #: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 City: \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_  
 Contact Name: \_\_\_\_\_  
 Phone: \_\_\_\_\_ Email: \_\_\_\_\_

### IMPRESSION RECOMMENDATIONS



- Use silicone impression material.
- Make impressions past the second bend.
- For optimal microphone placement, it is recommended to draw a horizontal line through the tragus on the earmold impression.

### PATIENT INFORMATION

Date: \_\_\_\_\_  
 First Name: \_\_\_\_\_  
 Last Name: \_\_\_\_\_

#### Audiometric Information

If no vent option is selected on this order form, Widex will select the appropriate vent based on the Audiometric Data provided.

	250Hz	500Hz (Required)	1kHz	2kHz	3kHz	4kHz
Right						
Left						

Widex may decrease the vent size to accommodate your order request. Would you like to be consulted before this change is made? (Y) (N)

Widex may change your selected sides for the volume control and program button features. Would you like to be consulted before this change is made? (Y) (N)

Previous User: (Y) (N) Age: \_\_\_\_\_  
 Ear Texture: ☐ Soft ☐ Medium ☐ Firm  
 Impression: ☐ Open mouth ☐ Closed mouth

☐ Use previous CAMISHA impression scan(s) to manufacture this order. Right Serial #: \_\_\_\_\_ Left Serial #: \_\_\_\_\_

### SELECT CUSTOM ITE OPTIONS

Model	Design	Hearing Instrument Colors	Venting May vary due to ear canal size	Additional Options	Canal Length
<input type="checkbox"/> EVOKE440 <input type="checkbox"/> EVOKE330 <input type="checkbox"/> EVOKE220 <input type="checkbox"/> EVOKE110 <small>ZEN Program comes standard in all EVOKE and UNIQUE models. (RC-Coil included. T-Coil not included. 312 battery standard.)</small>	Canal (R) (L) Half Shell (R) (L) Full Shell (R) (L)	<input type="checkbox"/> Beige <input type="checkbox"/> Medium Brown <input type="checkbox"/> Dark Brown <input type="checkbox"/> Clear <small>(If Clear shell is chosen, select a faceplate color)</small> <input type="checkbox"/> Beige <input type="checkbox"/> Medium Brown <input type="checkbox"/> Dark Brown	IROS Vent (R) (L) Semi IROS Vent (R) (L) Reverse Horn Vent (R) (L) Trench Vent (R) (L) Pressure Vent (R) (L) Largest Vent Possible (R) (L) No Vent (R) (L) XS (R) (L) S (R) (L) M (R) (L) L (R) (L) XL (R) (L)	Soft Hypoallergenic Shell Coat (R) (L) Hard Hypoallergenic Shell Coat (Default) (R) (L) Retention Ring (R) (L) Removal Line (R) (L) Removal Notch (R) (L) Canal Lock (N/A - Full and Half Shell) (R) (L) No Helix - 3/4 Shell (R) (L)	Medium (Standard) (R) (L) Short (R) (L) Long (R) (L) As Marked on Impression (R) (L)
Receiver	Feature*				Wax Guard
EU-IM (85 dB) (R) (L) EU-IP (90 dB) (R) (L)	Volume Control (R) (L) Program Button (R) (L) *To optimize size and ease of use, Widex recommends NOT adding a push button and volume control on both ears.				NanoCare™ (Standard) (R) (L) Extended Receiver Tubing (R) (L) Spring (R) (L)

### WARRANTY OPTIONS FOR ALL CUSTOMS (See EVOKE Price List for price information.)

INCLUDED WARRANTIES	OPTIONAL WARRANTIES
<b>EVOKE440/330</b> 3-Year Repair Warranty 3-Year Loss & Damage Coverage (Standard)	
<b>EVOKE220</b> 2-Year Repair Warranty 2-Year Loss & Damage Coverage (Standard)	3 <sup>rd</sup> Year Repair Warranty (R) (L) 3 <sup>rd</sup> Year of Loss & Damage Coverage (R) (L)
<b>EVOKE110</b> 1-Year Repair Warranty 1-Year Loss & Damage Coverage (Standard)	2 <sup>nd</sup> Year Repair Warranty (R) (L) 3 <sup>rd</sup> Year Repair Warranty (R) (L) 2 <sup>nd</sup> Year of Loss & Damage Coverage (R) (L) 3 <sup>rd</sup> Year of Loss & Damage Coverage (R) (L)

### SELECT ACCESSORIES (Indicate quantity)\*

TV-DEX: Silver _____	PHONE-DEX 2: Silver _____
ADD'L TV-DEX BASE: Silver _____	CALL-DEX: Black _____
COM-DEX: White _____	RC2-DEX: Silver _____
Grey _____	Black _____
Green _____	COM-DEX Remote Mic: Grey _____
UNI-DEX: Black _____	FM+DEX: Silver _____
	PerfectDry Lux Dryer: _____

\*All DEX devices are billable items for IM/IP models.

Additional Comments / Special Instructions: \_\_\_\_\_